

# Femoral Extraction Instruments

Designed to help remove various types of femoral implants



Loop  
Extractor



J-Hook  
Extractor



One-Piece  
Extractor



Standard Slap Hammer

PRODUCT NO'S:	
S1202	[Loop Extractor with Standard Slap Hammer]
S1202-01	[Loop Extractor Only] Overall Length: 6.5" (16,5 cm)
S1203	[J-Hook Stem Extractor with Standard Slap Hammer]
S1203-01	[J-Hook Stem Extractor Only] Overall Length: 4.75" (12,1 cm)
S1204	[One-Piece Stem Extractor with Standard Slap Hammer]
S1204-01	[One-Piece Stem Extractor Only] Overall Length: 4.125" (10,5 cm)
3925	[Standard Slap Hammer] 3/8"-16 Thread Gauge
3935	[Extra Large Slap Hammer] 3/8"-16 Thread Gauge



[www.innomed.net](http://www.innomed.net)

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362

© 2023 Innomed, Inc.