

# Femoral Extraction Instruments

Product No's:

S1202 [Loop Extractor with Standard Slap Hammer]

S1202-01 [Loop Extractor Only]

**Overall Length: 6.5"**

S1203 [J-Hook Femoral Stem Extractor with Standard Slap Hammer]

S1203-01 [J-Hook Femoral Stem Extractor Only]

**Overall Length: 4.75"**

S1204 [One-Piece Stem Extractor with Standard Slap Hammer]

S1204-01 [One-Piece Stem Extractor Only]


**Overall Length: 4.125"**

3925 [Standard Slap Hammer]

**Thread Gauge: 3/8"-16**

3935 [Extra Large Slap Hammer]

**Thread Gauge: 3/8"-16**



**FREE TRIAL**  
ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

**Phone** 800 548 2362  
912 236 0000  
**Fax** 912 236 7766  
**Internet** [www.innomed.net](http://www.innomed.net)  
**eMail** [info@innomed.net](mailto:info@innomed.net)

**Innomed**  
ORTHOPEDIC INSTRUMENTS  
103 Estus Drive Savannah GA 31404

