

Femoral Tibial Spreader with 90° Long Arms, Round Pads, and Measuring Scale

Designed by Adam Almaguer, MD


Longer 90° arms are helpful with visualization and useful for large patients

PRODUCT NO:

1866-D

Overall Length: 10" (25,4 cm)
Length to Bend: 9" (22,9 cm)
Depth from Bend: 2.75" (7 cm)
Pad Diameter: 1" (2,54 cm)

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY



Measuring scale in millimeters on the spreader handle

Reference Measurement Statement (Non-Measuring Function)

This device is not intended to provide quantitative measurements and does not incorporate a measuring function as defined under Regulation (EU) 2017/745 (MDR). Any scales, markings, or visual indicators present on the device are provided solely for reference and orientation purposes to support general use. These indicators are not calibrated against traceable standards and have not been validated for accuracy, precision, or repeatability.

Accordingly, values derived from these markings shall be considered approximate only and must not be used where accurate or traceable measurements are required for clinical decision-making, diagnosis, or other applications requiring defined metrological performance.

Users shall rely only on devices specifically intended and validated for measurement when quantitative accuracy is required.

www.innomed.net



FREE TRIAL ON MOST INSTRUMENTS

INNOMED 

ISO 13485:2016



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

© 2026 Innomed, Inc.