

Flared Darrach-type Bent Elevators

Designed by Richard L. Stowell, MD

Designed with a flare on the side to help further protect the deltoid

Left
#1967-L

Right
#1967-R

PRODUCT NO'S:	
1967-00	[Set]
Set Includes / Available Individually:	
1967-L	[Left]
Overall Length: 10.125" (25,7 cm)	
Blade Width: 1" 25,4 mm	
Blade Width with Flare: 1.25" (31,7 mm)	
Blade Depth: 3.5" (8,9 cm)	
1967-R	[Right]
Overall Length: 10.125" (25,7 cm)	
Blade Width: 1" 25,4 mm	
Blade Width with Flare: 1.25" (31,7 mm)	
Blade Depth: 3.5" (8,9 cm)	



www.innomed.net



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362