


Flexible Osteotome System

Provides an assortment of osteotome blades for various orthopedic surgery procedures

- ▶ Sharp, flexible blades are well suited for loosening implants from cement or bony ingrowth fixation
- ▶ Various blade widths and profiles allow great flexibility to follow the implant contours
- ▶ Modular handles are made of high impact surgical stainless steel and have a quick-coupling positive locking mechanism for ease of use and quick blade changes
- ▶ Slap hammer threads into the handle and is designed to facilitate blade removal

PRODUCT NO'S:	
S0011-00	[Complete Set]
Individual Instruments:	
S1002	[Thin Blade] 3" x 8mm
S1003	[Thin Blade] 3" x 10mm
S1004	[Thin Blade] 3" x 12mm
S1005	[Thin Blade] 3" x 20mm
S1006	[Curved Thin Blade] 3" x 12mm
S1007	[Curved Thin Blade] 3" x 20mm
S1008	[Thin Blade] 5" x 10mm
S1009	[Thin Blade] 5" x 8mm
S1020	[Handle with Quick-Coupling End] 6" <i>One Handle Only with this Product Number</i>
S1133	[Radial Blade] 5" x 10mm
S1120	[Radial Blade] 5" x 12mm (not shown)
S1134	[Radial Blade] 5" x 14mm
S1121	[Radial Blade] 5" x 16mm
S1122	[Radial Blade] 5" x 20mm (not shown)
S2007	[Small Slap Hammer] 12"
9018	[Case] (not shown)
Optional Blades (Not Included in Complete Set)	
S1123	[Extra Long Blade] 9" x 8mm
S1135	[Radial Blade Medial Curve] 6.75" x 11mm
S1136	[Radial Blade Lateral Curve] 6.75" x 11mm
S1137	[Radial Blade Medial Curve] 5" x 11mm
S1138	[Radial Blade Lateral Curve] 5" x 11mm



PRODUCT NO'S:	
S1222	[Chisel Blade] 2.5" x 8mm
S1223	[Chisel Blade] 2.5" x 10mm
S1224	[Chisel Blade] 2.5" x 12mm
S1225	[Chisel Blade] 2.5" x 20mm
S1228	[Chisel Blade] 5" x 10mm
S1229	[Chisel Blade] 5" x 8mm
S1230	[Chisel Blade] 5" x 20mm
S1231	[Chisel Blade] 5" x 12mm
S1020	[Handle with Quick-Coupling End] 6"

Flexible Chisel Blades

Sharp flexible blades are suitable for loosening implants from cement or bony ingrowth fixation

NEW



FREE TRIAL ON MOST INSTRUMENTS

© 2011 Innomed, Inc.

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

1.800.548.2362