

# Foster Cement Osteotome

Designed by Scott A. Foster, MD

*Designed to help remove UKA/TKA component*

Features a large handle and striking platform. The ultra hard titanium nitride coating helps to extend life by increasing surface hardness, prolonging sharpness, and resisting chemicals and corrosion.

**PRODUCT NO:**

5232

Osteotome Width: 6.7 mm

Overall Length: 8.5" (21,6 cm)

Handle Length: 5.75" (14,6 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
Fax +41 41 740 67 71

1.800.548.2362

# INNOMED