

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2005 Innomed, Inc., All Rights Reserved.

**Innomed**

103 Estus Drive, Savannah, GA 31404  
E-Mail: info@innomed.net



FREE TRIAL ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION

## Fromm Femur & Tibia Triangles

Designed to position and hold the femur and tibia during intramedullary nailing of the tibia, ligament repairs and extremity fractures. Allows knee to be flexed greater than 90° to allow reaming and nail insertion without displacing fracture. The triangles are available in three heights, 11", 14", and 16". The two smaller triangles are designed to fit inside the larger triangle for storage. They are supplied with an autoclavable silicone cushioning pad and velcro\* straps. The triangles are also radiolucent and gas or steam sterilizable.

Designed by  
Stuart E. Fromm, MD

Product No's: 2760-00 [Set of Three]  
2760-01 [11"]  
2760-02 [14"]  
2760-03 [16"]

\*Velcro® is a registered trademark  
of the Velcro Companies.

*Used for femur and tibia positioning during nailing, repairs and fractures*




— Over —




# Innomed




## Fromm Femur & Tibia Triangles




Tibial Nailing



Retrograde Femoral Nailing



Retrograde Femoral Nailing  
Triangle holds femur reduced  
(prevents sagging)



- Triangle Holds Tibia Reduced For:
- Open Reduction and Internal Fixation (ORIF) of tibia
  - Application of uniplane or multiplane external fixator
  - Knee ligament repairs and/or reconstruction