


Gap Clamp for Cortical Button Distal Bicep Repair

Designed by Corey Trease, MD

*Designed to be used to help consistently set
the gap for the radius cortex between the distal
biceps stump and the cortical button*



INNOMED

08/20 GH

PRODUCT NO:

5262

Overall Length: 8" (20.4 cm)

Jaw Width Tapers from: 6 to 2.5 mm

MADE EXCLUSIVELY
FOR INNOMED IN
GERMANY



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362

INNOMED