

Gelbke Freer Cement Trimmer/ Nerve Hook with TiN Coating

Designed by Martin K. Gelbke, MD

Designed to facilitate cement removal during total and partial knee replacement

- ▶ Consists of a freer elevator on one end and a nerve hook on the other
- ▶ Nerve hook accesses "tough to reach" corners of the knee
- ▶ Particularly useful for use with an ultra-congruent polyethylene insert, where trial liners are typically not used, once the final components have been placed
- ▶ Ultra hard titanium nitride coating helps to extend life by increasing surface hardness, prolonging sharpness, and resisting chemicals and corrosion



PRODUCT NO:

5007

Overall Length: 9.25" (23,5 cm)
Blade Width at End: 5 mm
Hook Depth: 5 mm

**PROUDLY
MADE
IN THE
USA**

ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED The Innomed logo features the word "INNOMED" in a bold, sans-serif font. To the right of the text is a stylized lowercase "i" inside a circle, with an American flag graphic integrated into the design.

1.800.548.2362