

Goldstein Glenoid Neck Retractor

Placed along the glenoid rim during open Bankart procedure to allow excellent exposure

The convex teeth sit easily into the glenoid rim while the strong end of the shaft allows the instrument to stay out of the surgeons view



PRODUCT NO:

T1004

Blade Width at Teeth: 18 mm

Blade Width at Widest: 36 mm

Overall Length: 8.5" (21,6 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

www.innomed.net



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

INNOMED

1.800.548.2362