


# Gupta Disc Space Spreaders with Easy Release Locking Mechanism

Designed by Munish C. Gupta, MD

*Designed to distract open collapsed disc spaces*

Locking ratchet mechanism helps prevent accidental release, and provides for controlled adjustment and easy release.



## PRODUCT NO'S:

4290 [Small Straight]  
Overall Length: 8.5" (21,6 cm)  
Blade Width: 10 mm  
Blade Thickness - Closed: 3 mm  
Opens to: 22 mm

4291 [Large Straight]  
Overall Length: 11.5" (29,2 cm)  
Blade Width: 13 mm  
Blade Thickness - Closed: 4 mm  
Opens to: 25 mm

4292 [Small Curved]  
Overall Length: 9.25" (23,5 cm)  
Blade Width: 10 mm  
Blade Thickness - Closed: 3 mm  
Opens to: 22 mm

4293 [Large Curved]  
Overall Length: 12.5" (31,8 cm)  
Blade Width: 13 mm  
Blade Thickness - Closed: 4 mm  
Opens to: 25 mm



[www.innomed.net](http://www.innomed.net)

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74  
info@innomed-europe.com

1.800.548.2362