

Hand/Finger Positioner

Designed by Emad Aboujaoude, MS, MPAS, PA-C

Designed to help provide surgical positioning during flouroscopy and fixation by isolating the operative digit while retracting the unaffected digits

Radiolucent positioner can be steam or gas sterilized.

Uses include but not limited to:

- ▶ Intramedullary Metacarpal Screw
- ▶ Phalanges CRPP
- ▶ Digit Amputation
- ▶ Digit Mass Excision
- ▶ Finger Joint Fusion




PRODUCT NO:

1134

Overall Length: 5.75" (14,6 cm)

Handle Width: 4.25" (10,8 cm)

Blade Width: 1.6" (4 cm)



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

INNOMED

1.800.548.2362