

# Hannum Rongeur/Grasper

Designed by Scott Hannum, MD

## *Teeth in jaw firmly holds bone and tissue*



Non-locking design can be easily gripped while allowing greater pressure to be applied. Available in three jaw sizes: short jaw for holding bone, medium jaw for smaller bones, and long jaw for tissue.

### PRODUCT NO'S:

1775-01 [Short Jaw]  
8mm Jaw Width  
Overall Length: 9.25"

1775-02 [Medium Jaw]  
5mm Jaw Width  
Overall Length: 9.25"

1775-03 [Long Jaw]  
3mm Jaw Width  
Overall Length: 9.25"



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax



# INNOMED

1.800.548.2362