

# Hohmann Retractor

*Designed like the original Hohmann-style retractor – made in the U.S.A.*

The OrthoLucent™ version is made of a lightweight carbon fiber material, which is completely radiolucent, helps to prevent from scratching component surfaces, and can be steam sterilized.



PRODUCT NO:
4558 [Standard] Overall Length: 9.5" Blade Width: 16mm
4558-R [OrthoLucent™] Overall Length: 9.625" Blade Width: 16mm

OrthoLucent™ is a trademark of MedX Composites

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73  
Fax +41 41 740 67 74

**INNOMED**

1.800.548.2362