



Modified Hohmann Retractors

Handle is contoured to allow better leverage and visualization

Useful for retracting tissues at the margins of the joint. Can be passed over the margins of the joint and held in place with weights or by hand.

The OrthoLucent™ version is made of a lightweight carbon fiber material, which is completely radiolucent, helps to prevent from scratching component surfaces, and can be steam sterilized.



PRODUCT NO'S:
6595 [Wide] Overall Length: 10" Blade Width: 41mm
4535 [Narrow] Overall Length: 10" Blade Width: 14mm
4535-R [OrthoLucent™ Narrow] Blade Width: 18mm Overall Length: 10"
4545 [Short-tipped Narrow] Designed by Carl DiRaimondo, MD Overall Length: 9.5" Blade Width: 14mm

OrthoLucent™ is a trademark of MedX Composites

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74



1.800.548.2362