

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2005 Innomed, Inc., All Rights Reserved.

Innomed

103 Estus Drive, Savannah, GA 31404
E-Mail: info@innomed.net



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Humeral Head Retractor

18cm length.

Product No:
T1007



Placed between the glenoid and the humeral head to obtain excellent exposure



Capsule Retractors

The convex teeth sit easily into the glenoid rim while the strong end of the shaft allows the instrument to stay out of the surgeons view.

25cm length.

Product No's:
T1008-01 [3 Prongs]
T1008 [2 Prongs]
T1009 [1 Prong]

Retractors for use in Bankart surgery

