

Humeral Head Retractor

Placed between the glenoid and the humeral head to obtain excellent exposure



PRODUCT NO:

T1007

Blade Width: 33 mm

Prong Width: 6 mm | 21 mm Gap | 6 mm

Overall Length: 7" (17,8 cm)

MADE FOR INNOMED IN
GERMANY

www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

INNOMED

1.800.548.2362