

Hur Modified Wide Mueller-type Femoral Neck Elevator

Designed by John Hur, MD

*Designed for the anterior approach, the wide design
helps to reduce stress on the proximal femur*



PRODUCT NO:

3416

Overall Length: 13" (33 cm)

Handle Length: 6.5" (16,5 cm)

Blade Width at Widest: 1.25" (31,7 mm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

www.innomed.net



INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362