

Incavo Wire Passer

Designed by Stephen J. Incavo, MD


Designed to pass multiple cerclage wires around a bone during a multiple wire wrap procedure



PRODUCT NO'S:

8610-01 [Small]
Overall Length: 7.5"
Accepts Wire Up To: 4mm (5/32")
8610-02 [Large]
Overall Length: 8.675"
Accepts Wire Up To: 4mm (5/32")

MADE
IN USA



FREE TRIAL ON MOST INSTRUMENTS

© 2011 Innomed, Inc.

INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

1.800.548.2362