

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2005 Innomed, Inc., All Rights Reserved.

Innomed

103 Estus Drive, Savannah, GA 31404
E-Mail: info@innomed.net



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Incavo Wire Passer

Designed to pass multiple cerclage wires around a bone during a multiple wire wrap procedure. Accepts a $5/32$ " diameter wire.

Designed by
S. Incavo, MD

Product No: 8610-02

Used for passing multiple cerclage wires around bone

