

# Joint, Calcaneal, Small Bone Compressor/Distractor

*Direction switches reverse the mechanism  
from compression to distraction*

Two hole sizes for pin size selection. Calibrated ratchet handle.


**PRODUCT NO:**

4865-LS

Overall Length: 8.5"

Holes For: .062" & .094" (1.6 & 2.5 mm) K-wire Pins

MADE  
IN USA



FREE TRIAL ON MOST INSTRUMENTS

© 2011 Innomed, Inc.

**INNOMED**



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 73  
Fax +41 41 740 67 74

1.800.548.2362