

Joint, Calcaneal and Small Bone Distractors

Two hole sizes and two arm designs allow for easier pin size selection and helps with distraction in a variety of indications



PRODUCT NO'S:	
Outspread Arms 4210-SB [Small] Hole Diameters: For 1/16" (1.6mm) & 1/8" (2.5mm) K-wire Pins Overall Length: 6"	Closed Arms 4210-SS [Small] Hole Diameters: For 1/16" (1.6mm) & 1/8" (2.5mm) K-wire Pins Overall Length: 6"
4210-LB [Large] Hole Diameters: For 1/16" (1.6mm) & 1/8" (2.5mm) K-wire Pins Overall Length: 8"	4210-LS [Large] Hole Diameters: For 1/16" (1.6mm) & 1/8" (2.5mm) K-wire Pins Overall Length: 8"

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
 www.innomed.net info@innomed.net

912.236.0000 Phone
 912.236.7766 Fax



INNOMED

1.800.548.2362