

Kaminsky OrthoLucent™ Browne-type Deltoid Retractors

Designed by Sean B. Kaminsky, MD

Used for the Delto-Pectoral approach, and can also remain in place for fracture reduction, plate positioning, and screw/wire/drill location confirmation

Used for acromioplasty, rotator cuff repair, and fracture fixation. Also helps to prevent from scratching component surfaces.

PRODUCT NO'S:

1670-01R [Small] Blade Width: 45mm Overall Length: 11.5"	1670-02R [Large] Blade Width: 57mm Overall Length: 11.5"
--	--

OrthoLucent™ is a trademark of MedX Composites



Completely radiolucent



Made of lightweight carbon fiber material—completely radiolucent and can be steam sterilized

FREE TRIAL ON MOST INSTRUMENTS

© 2011 Innomed, Inc.



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

1.800.548.2362