

Kaminsky Radiolucent Deltoid Retractors


Designed by Sean B. Kaminsky, MD

Designed for proximal humerus fractures—can remain in place for fracture reduction, plate positioning, and screw/wire/drill location confirmation

Helps to reduce operative time, assist in fracture reduction, and maintain hardware position without the frequent need for retractor removal and reintroduction.

PRODUCT NO'S:

1675-01 [Small] Opening Diameter: 22mm Overall Length: 11.75" Handle Length: 9.25" Blade Thickness: 5mm Blade Depth: 2.5"	1675-02 [Large] Opening Diameter: 28mm Overall Length: 11.75" Handle Length: 9.25" Blade Thickness: 5mm Blade Depth: 2.75"
--	---



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED

1.800.548.2362