



# Knee Retractors with Easy Grip Handles

*Helps provide excellent visibility and  
ligament protection during total and  
unicondylar knee replacement surgery*

Silicone handles help reduce holding fatigue.



Small Hohmann Retractor



Condylar Retractor



Soft Tissue Retractor



Superior Retractor

#### PRODUCT NO'S:

SS3035 [Small Hohmann Retractor]

Overall Length: 7"

Blade Width: 25mm

SS3037 [Condylar Retractor]

Overall Length: 7"

Blade Width: 12mm

SS3038 [Soft Tissue Retractor]

Overall Length: 8.25"

Blade Width: 36mm

SS3042 [Superior Retractor]

Overall Length: 8.25"

Blade Width: 31mm



FREE TRIAL ON MOST INSTRUMENTS

© 2011 Innomed, Inc.

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 73  
Fax +41 41 740 67 74

1.800.548.2362