

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2005 Innomed, Inc., All Rights Reserved.

Innomed

103 Estus Drive, Savannah, GA 31404
E-Mail: info@innomed.net



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Kolbel Self-Retaining Glenoid Retractor with Center Blade

Two pairs of snap-in, freely pivoting blades included.

Product No's:
T1050 [Complete Set]
T1050-01 [Retractor Only]
T1050-02 [Center Blade]
T1018-P [36mm X 36mm Blades - Pair]
T1019-P [36mm X 53mm Blades - Pair]



Blade can be reversed for shallow or deep retraction



Kolbel Self-Retaining Retractor Blades

Product No's:
T1018 [36mm X 36mm Blade]
T1019 [36mm X 53mm Blade]
T1020 [36mm X 68mm Blade]
T1021 [36mm X 85mm Blade]

Product No's:
T1022 [20mm X 36mm Blade]
T1023 [20mm X 53mm Blade]
T1024 [20mm X 68mm Blade]
T1025 [20mm X 85mm Blade]

