

1-800-548-2362

International (912) 236-0000




www.innomed.net

(912) 236-7766 Fax

©2005 Innomed, Inc., All Rights Reserved.

# Innomed

103 Estus Drive, Savannah, GA 31404  
E-Mail: info@innomed.net



FREE TRIAL ON MOST INSTRUMENTS  
CALL FOR MORE INFORMATION

## Lambotte Osteotomes with Handle

Designed with a handle for better control, which helps reduce rotation of the osteotome during use. The handle also provides a larger striking area for use with a mallet. The osteotome shafts are manufactured with stainless steel and are available both straight and curved.

Designed by  
J. Cherf, MD

Product No's: 5250-01 [1/4" Straight]  
5260-01 [1/4" Curved]

*Handle allows for better control,  
reducing rotation during use*

