

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2000 Innomed, Inc., All Rights Reserved.

Innomed
Instrument Division
103 Estus Drive, Savannah, GA 31404



FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Large Bone Curettes

Designed with a 5/16" diameter shaft allowing better visualization into the medullary canal during use. The contoured handle is designed to keep the curette from slipping in the surgeon's hand and for better control. Shaft length 10-1/2".


The Angled Large Curette is designed for use in the acetabulum or exposed bone. The shaft diameter is 5/16" and has a contoured handle. Shaft length 10-1/2".

Product No's: 5160 [Set of 1 each including Case]
5160-01 [Angled Small]
5160-02 [Straight Small]
5160-03 [Angled Medium]
5160-04 [Angled Large]
5160-05 [Straight Medium]
Storage Case Available

Allow better visualization into the medullary canal



Storage Case Available



Curette Ends at Actual Size

