

Lawton Flexor Tendon Repair Clamp




Helps the surgeon to reapproximate the lacerated flexor tendon and tie the core suture(s)

Luer lock ends allow the interchange of two 25g disposable hypodermic needles.

PRODUCT NO:

1795
Overall Length (flat): 6"
Overall Length (bent): 4.5"
Leg Length (without needle): 1.25"

Designed by Jeffrey Lawton, MD



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED

1.800.548.2362