

# Lighted Hip Retractors

*Lighting attachment for enhanced visual exposure*

Excellent for use in acetabular exposure and total hip replacements.

Especially useful for anterior approach.

Can be attached to a fiber optic light cable with ACMI (female) connector. Retractors can be steam sterilized.



*Jana Lighted  
Cobra Retractors*

Designed by Ajoy K. Jana, MD

*Lighted Standard  
Cobra Retractor*

#### PRODUCT NO'S:

**6119-L** [Jana Lighted Cobra]

Overall Length: 14.75" (37,5 cm)

Blade at Widest: 33 mm

**6130-L** [Lighted Standard Cobra]

Overall Length: 12" (30,5 cm)

Handle Length: 7" (17,8 cm)

Blade Width at Widest: 33 mm



## Light Source Cable Adapters

#### PRODUCT NO'S:

**8009-S** [ACMI to Storz Adapter]

**8009-W** [ACMI to Wolf Adapter]



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

[www.innomed.net](http://www.innomed.net)



# INNOMED



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**