




Now Available with
Diamond Cut Pads


Lombardi Gap Balancing Femoral Tibial Spreader

Designed by Adolph V. Lombardi Jr., MD

Designed to help separate the femur and tibia during total knee procedures, with the pads being parallel when measured at 20mm of separation



PRODUCT NO'S:	
Horizontal Grooved Pads	Diamond Cut Pads
1878 [Large] Overall Length: 9.25" Pads: 22mm x 13mm Opens to 50mm	1878-D [Large] Overall Length: 9.25" Pads: 22mm x 13mm Opens to 50mm
1877 [Small] Overall Length: 7" Pads: 22mm x 13mm Opens to 35mm	1877-D [Small] Overall Length: 7" Pads: 22mm x 13mm Opens to 35mm



FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED

1.800.548.2362