

Long Handled Fukuda Type Retractors

Design modification of an original design by Evan Flatow, MD & Louis Bigliani, MD

Designed to retract the humeral shaft posteriorly, helping to expose the entire glenoid surface, these retractors feature a longer handle for use when added leverage and depth is needed for use in large patients



Narrow with Long Handle
#1930-H

Wide with Long Handle
#1940-H

New!
Larger
Diameter
Handle

PRODUCT NO'S:

1930-H [Narrow with Long Handle]

Blade Width: 32 mm

Opening: 25 x 40 mm

Overall Length: 9.35" (23,7 cm)

1940-H [Wide with Long Handle]

Blade Width: 38 mm

Opening: 32 x 40 mm

Overall Length: 9.35" (23,7 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



www.innomed.net

INNOMED



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362