


Lower Extremity Leg Positioner

Designed by Ronald Romanelli, MD

Designed to lift the knee for lower extremity casting applications



Supplied with one autoclavable silicone pad.
Positioner is radiolucent and gas or steam sterilizable.

PRODUCT NO'S:

2745

Dimensions: 5.5" H x 9.5" L x 9.25" W

Replacement Parts:

2760-P [Silicone Pad]

MADE
IN USA

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

© 2011 Innomed, Inc.

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

INNOMED

1.800.548.2362