


# MIS Utility Knee Retractor

Designed by William Robb, MD

*Used interchangeably for medial exposure, lateral exposure  
and to assist in posterior exposure for the tibia*

Helps to keep hands out of the field of view while providing retraction in minimally invasive knee surgery.



**PRODUCT NO:**

3220-03

Overall Length: 9" (22,9 cm)

Blade Width: 1,6 cm



ISO 13485:2016

**FREE TRIAL ON MOST INSTRUMENTS**

[www.innomed.net](http://www.innomed.net)

**INNOMED** 



103 Estus Drive, Savannah, GA 31404  
info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax

**Innomed-Europe** Tel. +41 41 740 67 74  
info@innomed-europe.com

**1.800.548.2362**