


| PRODUCT NO'S: | |
|-------------------------|--|
| S7500-01 | [Complete Set with Case] |
| Individual Instruments: | |
| S7505 | [Narrow Cement Removal Gouge, Short] Shaft Length: 15 cm Gouge: 9 mm, negative |
| S7507 | [Narrow Cement Removal Gouge, Long] Shaft Length: 24 cm Gouge: 9 mm, negative |
| S7510 | [Narrow Offset Cement Removal Gouge] Shaft Length: 24 cm Gouge: 9 mm, negative |
| S7515 | [Acetabular Chisel] Shaft Length: 24 cm Chisel: 7.5 mm |
| S7520 | [Offset Chisel] Shaft Length: 15 cm Chisel: 9 mm |
| S7525 | [Flared Angle Gouge] Shaft Length: 24 cm Gouge: 9 mm, positive, angle 15° down |
| S7530 | [Wide Gouge] Shaft Length: 24 cm Gouge: 11.5 mm, negative |
| S7535 | ["V" Splitter] V-Shaped Chisel: 7 mm |
| S7587 | [Saddle Punch] Shaft Length: 24cm Punch: 16.5 mm x 6.5 mm |
| S7590 | [Cement Splitting Osteotome] Shaft Length: 24 cm |
| S7595 | [Cement Removal Osteotome, Short] Shaft Length: 15 cm Osteotome: 8 mm |
| S7597 | [Cement Removal Osteotome, Long] Shaft Length: 24 cm Osteotome: 8 mm |
| S7540 | [4.4 mm Drill] |
| S7545 | [4.4 mm Drill Guide] |
| S7550 | [6.4 mm Drill] |
| S7555 | [6.4 mm Drill Guide] |
| S7560 | [Straight Cement Removal Hook] Hook Curette: 10 mm |
| S7565 | [Curved Cement Removal Hook] Hook Curette: 10 mm |
| S7570 | [Cross Bar] |
| S7575 | [7 mm T-Handle Conical Tap] |
| S7580 | [9 mm T-Handle Conical Tap] |
| S7585 | [Slotted Mallet] |
| 9075 | [Case Only] |




Mueller-Type Cement Removal Instruments



Used for cement removal in the knee, hip, and shoulder


Optional Instruments:



T-Handle Chuck for use with Drills

| PRODUCT NO'S: | |
|---------------|------------------------|
| 8247-00 | [T-Handle Chuck & Key] |
| 8247-01 | [T-Handle Chuck] |
| 8247-02 | [Chuck Key] |

Conical Tap & Mallet



Drill & Drill Guide



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

1.800.548.2362