


Mobasser Bayonet Forcep

Designed so that the tips can hold a larger amount of bone graft, and then be used as a tamp to push it down into the disc space

Bayonetted forcep has a long reach, which is particularly helpful in less invasive spine surgery.



Designed by
Jean Pierre Mobasser, MD

Product No: 5082
OVERALL LENGTH: 10"
FORCEP ENDS: 5x10MM

FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2008 Innomed, Inc., All Rights Reserved

Phone 800 548 2362
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

