

Modified Closed Blade Fukuda Retractor

Designed by Clint Merritt, PA-C

Designed to help position the humerus posterior and inferior to the glenoid during open shoulder reconstruction or shoulder replacement procedures



PRODUCT NO:
1945
Overall Length: 8.25" (21 cm)
Blade Width: 1.25" (3,2 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

www.innomed.net



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

INNOMED

1.800.548.2362