



Modified Lambotte Cup Removal Osteotomes

Designed with different hemisphere of curves to match cups of different sizes

Four osteotomes with different hemispherical radii allow the osteotomes to fit next to the outer surface of different size acetabular hip cups. The osteotomes have a handle for better control and a hammering platform.



| PRODUCT NO'S: |
|---|
| 5240-44 Blade Radius: 44mm Overall Length: 12.75" Handle Length: 4.75" |
| 5240-48 Blade Radius: 48mm Overall Length: 12.75" Handle Length: 4.75" |
| 5240-52 Blade Radius: 52mm Overall Length: 12.75" Handle Length: 4.75" |
| 5240-56 Blade Radius: 56mm Overall Length: 12.75" Handle Length: 4.75" |

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax



INNOMED

1.800.548.2362