


Modified Fukuda-type Retractor with Reamer Slot

Designed by Richard J. Miller, MD

Center cutout slot allows the shaft of a reamer to fit more posteriorly

Used to retract the humeral shaft posteriorly and helping to expose the entire glenoid surface. Available in small and large sizes.



PRODUCT NO'S:

1898 [Narrow]

Overall Length: 7.25" (18,4 cm)

Blade Width: 32 mm

Opening: 25 x 40 mm

1899 [Wide]

Overall Length: 7.25" (18,4 cm)

Blade Width: 38 mm

Opening: 32 x 40 mm



www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362

INNOMED