





Left
#3012-L



Center
#3012-01

Not Included In Set.
Available Separately.



Right
#3012-R

Modified Wide Hohmann Retractor with Taylor Tip

Designed by Jeffrey P. Beckenbaugh, DO

Anterior and posterior acetabular retractors for all approaches, including the direct anterior approach, featuring a hammer platform for insertion with a mallet

Used as a calcar and posterior femoral retractor for the posterior approach, and an anterior femoral elevator for the direct anterior approach.



Left
#3012-L

Right
#3012-R

Center
#3012-01

Not Included In Set.
Available Separately.

PRODUCT NO.'S:	
3012-00	[Set - Left & Right]
Set Includes / Available Individually:	
3012-L	[Left]
Overall Length: 14.75" (37,5 cm)	
Blade Depth: 3.25" (8,3 cm)	
Blade Width: 45 mm	
3012-R	[Right]
Overall Length: 14.75" (37,5 cm)	
Blade Depth: 3.25" (8,3 cm)	
Blade Width: 45 mm	
Not In Set / Available Separately:	
3012-01	[Center]
Overall Length: 14.75" (37,5 cm)	
Blade Depth: 3.25" (8,3 cm)	
Blade Width: 45 mm	



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362



www.innomed.net

INNOMED