

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2005 Innomed, Inc., All Rights Reserved.

Innomed

103 Estus Drive, Savannah, GA 31404
E-Mail: info@innomed.net




FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION

Modular Weights

Designed to hook onto a retractor handle which helps the retractor stay in place by the weights attached rather than being held. The desired weight can be adjusted by attaching additional weights to the primary weight which has the attaching hook.

Product No's:
3430-01 [1.5 lbs.]
3430-02 [2.0 lbs.]
3430-03 [2.5 lbs. with attaching hook]



Helps free assisting personnel

