

Non-Metallic Mallet

Designed by Justin Duke, MD

Designed to be used when adjustments are needed for implants. The non-metallic material insures that no metal is transferred to the implant.



PRODUCT NO:

7834

Overall Length: 8.25" (21 cm)

Handle Length: 4.75" (12 cm)

Head Diameter: 1.5" (3,8 cm)



ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

www.innomed.net



INNOMED 



103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362