


O'Brien Bone Clamp with Drill Guide

Designed by
Todd O'Brien, DPM

Calibrated handle measures bones from 6mm up to 24mm.
Guide accommodates 1.5, 2.0 and 2.7mm screws.

6.2mm diameter guide, 5" overall length

Product No: 1815



FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2007 Innomed, Inc., All Rights Reserved

Phone 800 548 2362
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

