

## O'Brien Nail Stat

*Combining the advantages of a hemostat and a spatula, this instrument facilitates atraumatic partial and total nail avulsion*

When used for partial avulsion, instrument may be inserted beneath the nail plate prior to cutting with a Beaver blade. The lower jaw of the clamp protects the nail bed from laceration as the blade is advanced.



**PRODUCT NO:**  
1565  
Overall Length: 6"  
Length from Bend : 3"

Designed by Todd O'Brien, DPM

FREE TRIAL ON MOST INSTRUMENTS



103 Estus Drive, Savannah, GA 31404  
www.innomed.net info@innomed.net

912.236.0000 Phone  
912.236.7766 Fax



# INNOMED

1.800.548.2362