

OrthoLucent™ Modified Fukuda-type Retractors

Used to retract the humeral shaft posteriorly, helping to expose the entire glenoid surface, the carbon fiber PEEK composite material is strong, lightweight, completely radiolucent, helps to prevent from marring component surfaces, and can be steam sterilized



OrthoLucent™ Wide
1940-R

OrthoLucent™ Narrow
1930-R



PRODUCT NO'S:

1930-R [OrthoLucent™ Narrow]
Blade Width: 32 mm
Opening: 25 x 40 mm
Overall Length: 7.5" (19,1 cm)

1940-R [OrthoLucent™ Wide]
Blade Width: 38 mm
Opening: 32 x 40 mm
Overall Length: 7.5" (19,1 cm)

MADE EXCLUSIVELY
FOR INNOMED IN
SWITZERLAND

www.innomed.net

ISO 13485:2016

FREE TRIAL ON MOST INSTRUMENTS

INNOMED

USA | CE | © 2023 Innomed, Inc.

103 Estus Drive, Savannah, GA 31404
info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
info@innomed-europe.com

1.800.548.2362