

OrthoLucent™ O'Brien Bone Clamp

Designed by Todd O'Brien, DPM

Designed for use in stabilization of a fracture

The carbon fiber PEEK material is strong, lightweight, completely radiolucent, can be steam sterilized, and helps to prevent from marring component surfaces.

PRODUCT NO:
1815-R
Overall Length: 5.25" (13,3 cm)

MADE EXCLUSIVELY
FOR INNOMED IN
SWITZERLAND



ISO 9001:2008 • ISO 13485:2003

FREE TRIAL ON MOST INSTRUMENTS

Scan to
Launch Our
Website



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 74
Fax +41 41 740 67 71

INNOMED

1.800.548.2362