

Offset Cup Liner Inserter

Offset to improve visualization and for mis hip surgery

Product No's:

5022 [22mm head dia.]
5026 [26mm head dia.]
5028 [28mm head dia.]
5032 [32mm head dia.]
5036 [36mm head dia.]
5038 [38mm head dia.]

*Overall Length: 16.25"
Shaft Diameter: .5" (12.7mm)*



FREE TRIAL
ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



©2009 Innomed, Inc., All Rights Reserved

Phone 800 548 2362
912 236 0000
Fax 912 236 7766
Internet www.innomed.net
eMail info@innomed.net

Innomed
ORTHOPEDIC INSTRUMENTS
103 Estus Drive Savannah GA 31404

