

1-800-548-2362

International (912) 236-0000



www.innomed.net

(912) 236-7766 Fax

©2002 Innomed, Inc., All Rights Reserved.

Innomed
Instrument Division
103 Estus Drive, Savannah, GA 31404

FREE TRIAL ON MOST INSTRUMENTS
CALL FOR MORE INFORMATION



Offset Punches for Total Hip Stem Removal

Two sizes of offsets allow the punches to be used to tap on a distal portion of an intramedullary nail, after a window has been made in the bone below the tip of the nail.

Product No's: 5125-01 [1/2" offset]
5125-02 [1-1/4" offset]

Used to help remove an intramedullary nail or prosthesis stem via a window in the shaft of the bone

