


Offset Punches

Helps in the removal of hip stems

Used to help remove a hip prosthesis stem via a window in the shaft of the femur. Two sizes of offsets allow the punches to be used to tap on a distal portion of the hip stem, after a window has been made in the femur below the tip of the stem.



PRODUCT NO'S:

5125-01 [1/2" Offset]
Overall Length: 11" (27,9 cm)
Punch End Offset: .5" (1,3 cm)
Punch End Diameter: 7 mm

5125-02 [1.25" Offset]
Overall Length: 11" (27,9 cm)
Punch End Offset: 1.25" (3,2 cm)
Punch End Diameter: 7 mm

