


OrthoLucent™ Hohmann Retractors

Carbon fiber material is lightweight, completely radiolucent and can be steam sterilized

Also helps to prevent from scratching component surfaces.

PRODUCT NO'S:
4535-R [Modified] Blade Width: 18mm Overall Length: 10"
4550-R [Modified Blunt] Blade Width at Widest: 24.5mm Overall Length: 10.75"
4558-R [Standard] Blade Width: 16mm Overall Length: 9.625"

Ortho Lucent™ is a trademark of MedX Composites



Modified Hohmann

Modified Blunt Hohmann

Standard Hohmann

FREE TRIAL ON MOST INSTRUMENTS

© 2011 Innomed, Inc.

INNOMED 



103 Estus Drive, Savannah, GA 31404
www.innomed.net info@innomed.net

912.236.0000 Phone
912.236.7766 Fax

Innomed-Europe Tel. +41 41 740 67 73
Fax +41 41 740 67 74

1.800.548.2362